ComStation^{CT} Flamenco, Best.-Nr. 77 0606 00

ComStation^{CT} Flamenco, order no. 77 0606 00

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- 1 - ComStation^{CT} Flamenco
- 2 - Connection plug
- 3 - * Connection socket ComStation

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

- Press the left and right function keys simultaneously (3 sec.) until the configuration menu is shown in the display.

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●	Circle key	Set the marked menu item.
▲	Upwards arrow key	Scroll upwards through the list.
▼	Downwards arrow key	Scroll downwards through the list.
⌂	Home key	Cancel, without storing.

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Menu items

Language:	Selecting the user language.
RAN > RAN number:	Setting of number of room devices (= RAN number).
RAN > Test RAN:	Test if room devices are ready to operate and if they are correctly connected to the RAN.
RAN > Status:	Check whether the stored RAN number is equal to the number of operational room devices.
Address:	Setting of room address.
Cleaning time:	Setting the time (seconds) how long it takes to wipe the front of the ComStation ^{CT} .
Contrast:	Function use only for Tunstall technicians.
Key sound:	Switching the key sound on or off. "Key sound ON" (factory setting) is recommended.
Fault tone:	Setting the fault tone volume.
Info:	Displaying the ComStation's software revision.
Audio Test:	Function use only for Tunstall technicians.
Key Function:	Setting the function of Alarm Key, Presence Key 1 and Presence Key 2 on the ComStation ^{CT} .
Reset:	Function use only for Tunstall technicians. (Restart the ComStation ^{CT}).

2. Entering of necessary settings**Selecting the user language**

1. Mark "Language" using the arrow keys. Then press the circle key.
2. Mark the desired language using the arrow keys: D = German; GB = English, etc.
3. Press the circle key to set the selection.

Setting of RAN number (0 – 30) (Room Area Network)

RAN number = Number of devices in the room (switches, connection sockets, room lamps, etc.) that are connected via RAN. (**No** pear push switch, **no** patient handset).

1. Mark "RAN" using the arrow keys. Then press the circle key.
2. Mark "RAN number" using the arrow keys. Then press the circle key.
3. Mark the desired RAN number using the arrow keys.
4. Press the circle key to set the selection.

Setting of room address (0 – 110)

1. Mark "Address" using the arrow keys. Then press the circle key.
2. Mark the desired room address using the arrow keys.
3. Press the circle key to set the selection.

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3. Entering of optional settings

The remaining menu items provide optional settings. Set these if required. Examples:

Fault tone volume soft, medium, loud or off

WARNING! The tone that attracts attention for a fault may only be turned off, if it is made sure, that the attention is attracted in another way.

1. Mark "Fault Tone" using the arrow keys. Then press the circle key.
2. Mark "Soft", "Medium", "Loud" or "Off" using the arrow keys.
3. Press the circle key to set the selection.

Setting the key function

CAUTION! Switched off keys are not available for nursing staff!

Alarm Key

With the factory settings the alarm key on the ComStation^{CT} is operational and is initiated by pressing it.

If there is the risk, that the alarm button is initiated inadvertently, it might make sense to set a delay time (2 or 3 seconds). In that case the alarm will be initiated after the alarm key has been pressed and then hold for 2 or 3 seconds.

CAUTION! Setting a delay time does not comply with the German standard DIN VDE 0834. A delay time may only be set, where compliance with the DIN VDE 0834 is not required.

In well-founded exceptional cases you can switch of the alarm key function.

NOTE! Die Alarm Key is only active while the staff presence is switched on in the room.

Presence Key 1 / Presence Key 2

With the factory settings the presence keys are operational. In exceptional cases it may be useful to switch of the key function for these keys.

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4. Checking the room bus RAN

Status

Check whether the set RAN number (refer to section “Setting of RAN number”) equals to the number of operational room devices.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Mark “Status” using the arrow keys. Then press the circle key.

Display: **OK**: The set RAN number is equal to the number of operational room devices.

Display: **Error**: The set RAN number is not equal to the number of operational room devices.

3. Press the home key to end. In case of an **error** perform a “Test RAN”.

Test RAN

Test if room devices are ready to operate and if they are correctly connected to the RAN.

Check the correct setting of the RAN number.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Test start: Mark “Test RAN” using the arrow keys. Then press the circle key.

The LEDs of all connected room devices must flash. Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.

3. Compare the number of flashing room devices with the set RAN number. If the set RAN number is wrong, correct the RAN number setting after the RAN test as described in section “Setting of RAN number”.
4. End of test: Press the home key.

5. Exit the configuration menu

When all settings are made and all tests are completed, you have to exit the configuration menu:

- Press home key several times until normal operation display appears.

NOTE! If for a period of three minutes no function key is pressed, the system will leave the configuration menu.

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ComStation^{CT} Flamenco, order no. 77 0606 00

Terminal with two-way speech communication, designed for desktop installation in the nurse station, incl. red Call Key, blue Alarm Key, green Presence Key (staff 1), yellow Presence Key (staff 2), 4 Function Keys as soft keys, microphone and loudspeakers.

NOTE! The complete installation of the nurse call system is described in the Technical Manual.


Glass – fragile! The operating front cover of the ComStation^{CT} is made of glass. Handle with care!

Installation

1. Connect the connection plug [2] of the ComStation^{CT} Flamenco [1] to the connection socket ComStation [3].
2. Tighten the two screws of the connection plug [2], to protect the connection plug against unintentional disconnection.

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- 1 - ComStation^{CT} Flamenco 3 - * Connection socket ComStation
2 - Connection plug

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

* Not included with ComStation^{CT} delivery.

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Symbols in the configuration menu

●	Circle key	Set the marked menu item.
▲	Upwards arrow key	Scroll upwards through the list.
▼	Downwards arrow key	Scroll downwards through the list.
🏠	Home key	Cancel, without storing.

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4. Checking the room bus RAN

Status

Check whether the set RAN number (refer to section "Setting of RAN number") equals to the number of operational room devices.

1. Mark "RAN" using the arrow keys. Then press the circle key.
2. Mark "Status" using the arrow keys. Then press the circle key.

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3. Press the home key to end. In case of an **error** perform a "Test RAN".

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3. Compare the number of flashing room devices with the set RAN number. If the set RAN number is wrong, correct the RAN number setting after the RAN test as described in section "Setting of RAN number".
4. End of test: Press the home key.

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