**IP-SystemManager, Bestell-Nr. 76 2100 00**

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**IP-SystemManager, order no. 76 2100 00**

The control of a Flamenco<sup>IP+</sup> or Flamenco<sup>IP</sup> nurse call system is performed by decentralised LAN networked IP-SystemManagers.

All IP-SystemManagers are equipped with identical hardware. However, each IP-SystemManager has a defined function. For this function the required software and a project specific database were factory installed. That means, each IP-SystemManager must be installed at that location, for which it has been prepared at the factory.

For installation location and IP address refer to the device label.

**NOTE!** The complete installation of the system is described in the Technical Manual.

**CAUTION!** The printed circuit board behind the service cover includes electrostatic sensitive components. Avoid touching.

**A Product description**

- 1 - Housing \* Not included in the scope of delivery, please order separately.
- 2 - Service cover
- 3 - \* DIN rail
- 4 - 2 DIN rail clamps

**B Mounting**

Installation in dry rooms only. Always easily accessible to authorised persons (access space at least 60 cm wide).


The IP-SystemManager must be protected against unintentional contact, e.g. by installing it in a 19" wall mounting cabinet, 7HU, GD (order no. 77 3519 07).

In case of installation in switch cabinets or similar, dissipated heat shall be discharged by means of forced ventilation, if necessary.

Mounting on DIN rail according to fig. B. e.g. with 19" mounting set (order no. 76 0900 01).

**C Dismantling**

Dismantling from the DIN rail according to fig. C.

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Patchkabel, mindestens Kat. 5e.

## Connections

★ The connections marked with a star are only enabled, if the IP-SystemManager was prepared accordingly at the factory, refer to the project specific documents.

### Connection "Power"

Connection to the 24 V power supply ring of the ward or directly to a 24 V power supply unit. Insert insertion bridges!

Conductor cross-section: 2.5 mm<sup>2</sup>

Connector, 4-pole: spare part order no. 00 0211 33.

Two insertion bridges: spare part order no. 00 0220 52.

### Connection "UPS Input"

Input for the signalling contact of a connected uninterruptible power supply (UPS), which signals buffer operation. This signalling contact is available at:

- UPS Control unit 10A (77 3411 00): "Buffering" contact.
- DC-UPC with capacitor storage 24V, 6kWs (77 3413 00): "Buffering" contact.
- Power supply unit UPS (77 3400 00): "Battery active" contact.
- Power supply unit UPS 60 (77 3400 10): "Battery active" contact.

Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>

Connector, 2-pole: spare part order no. 00 0211 47.

### ★Connection "Fault Relay"

The signalling of the output is in parallel with the "Fault" LED (below "Fault Relay"), for details refer to the chapter "LED displays".

Fault message relay output.

Change over contact, potential free.

Maximum contact load: 2 A, 24 V

Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>

Connector, 3-pole: spare part order no. 00 0211 45.

### ★Connections "OSYnet 1" and "NF 1"

Connection of a group bus OSYnet. For the maximum number of bus users that can be connected, please refer to the Planning manual.

**CAUTION!** The IP-SystemManager must not be installed in the middle of the group bus line. It must be connected at one end of the group bus line, as it includes an integrated bus terminating resistor.

In nurse call systems without speech communication, connection "NF 1" is not used.

For permitted cable types, refer to the Technical manual.

Connector, 3-pole, for OSYnet: spare part order no.

00 0211 45; Connector, 4-pole, for NF: Spare part order no. 00 0211 36.

### ★Connections "OSYnet 2" and "NF 2"

Not used.

### Connection "Ethernet"

Connection to the Ethernet of the nurse call system. This connects the IP-SystemManager with the other IP-SystemManagers and the IP based devices of the nurse call system.

**NOTE!** The Ethernet port is internally equipped with a separation point according to EN 60601-1 (2 x MOPP). An external separation device is not required.

Ethernet, 10/100 Mbit LAN connection, RJ45 socket.

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**\*Connection "RS232 A"**

ESPA 4.4.4 protocol.

For example, connection of a radio paging system, DECT system, fire alarm system or medical electrical system.

**NOTE!** If room terminals are connected to the OSYnet port and an RS232 port is used a separator device according to EN 60601-1 (2 x MOPP) must be connected between the connected device and the RS232 port, e.g. the interface isolator RS232 (order no. 76 5010 00).

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
**EN Connecting**

1. Strip the connection cables to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip off the wires 7 mm.

**CAUTION!** If the IP-SystemManager is operated without insertion bridges at the power connections, high operating currents from the system flow through the IP-SystemManager. This may damage the device. The two insertion bridges **must** be inserted according to fig. E!

2. Connect the wires to the connectors (pluggable screw-type connectors) according to fig. D.
3. Insert the patch cable's RJ45 plug into the Ethernet socket.

**LED displays****Green LED "Power"**

LED lights up permanently. = Power supply available at the IP-SystemManager.

**Green LED "Status"**

LED flashes once a second. = IP-SystemManager is ok.  
LED lights up permanently or is off. = Fault! IP-SystemManager is not operational.

**Red LED "Fault" below "Fault Relay"**

LED lights up permanently. = There is at least one fault. Signalling in parallel to "Fault Relay" output.

This is a collective indicator. The assigned fault locations depend on the operating mode (see device label) of the IP-SystemManager:

**- Operating mode "System" or "System + Local"**

Collective indicator for faults in the complete nurse call system. Indicated fault types:

- An OSYnet bus user is defective or disconnected.
- A RAN bus user is defective or disconnected.
- The power supply on a ward is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.
- The power supply at an IP-SystemManager is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.
- An OSYnet bus line is interrupted.
- LAN connection to an IP-SystemManager is interrupted.
- An IP-SystemManager is not operational.
- The power supply unit is in buffer mode. Indication in parallel to the red LED "Fault" below "UPS Input".

**- Operating mode "Local"**

Collective indicator for faults in the ward connected to this IP-SystemManager. Indicated fault types:

- An OSYnet bus user is defective or disconnected.
- A RAN bus user is defective or disconnected.
- The power supply unit is in buffer mode. Indication in parallel to the red LED "Fault" below "UPS Input".

**NOTE!** The red LED "Fault" (below "Fault Relay") and the "Fault Relay" connection of IP-SystemManagers with operating mode "BMA/MED", "Voice Gateway" and "Voice Gateway Master" are inactive.

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**Red LED "Fault" below "UPS Input"**

LED lights up permanently. = The uninterruptible power supply operates in buffer mode.

**OSYnet 1: Green LED "Rx", yellow LED "Tx"**

LEDs are flashing. = Data traffic on the group bus, connected to "OSYnet 1".

Yellow LED lights up permanently. = Group bus OSYnet 1 has not yet been configured with the SystemOrganizer.

**OSYnet 2: Green LED "Rx", yellow LED "Tx"**

Yellow LED lights up permanently. = There is no active use licence for OSYnet 2.

**Yellow LED "NF 1"**

LED lights up permanently. = There is a speech connection or an announcement in the group connected to OSYnet 1.

**Yellow LED "NF 2"**

Not used.

**Ethernet, green LED "Link", yellow LED "Speed"**

The LEDs are signalling "Link" and "Speed" of the Ethernet connection.

**RS232 A, green LED "Rx", yellow LED "Tx"**

LEDs are flashing. = Data traffic between the IP-SystemManager and the device connected to „RS232 A“.

**RS232 B, green LED "Rx", yellow LED "Tx"**


As for "RS232 A", but for RS232 B.

**F Service cover**

The service cover [2] is removed by levering off. After use re-attach the service cover.

- 5 - 6 grommets for fixing the service cover
- 6 - Li battery, 3 V, CR2032
- 7 - Reset button
- 8 - DIP switch for addressing: Do not change the factory setting!

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**Replace battery**

Behind the service cover [2] there is a lithium button cell [6] (3 V, CR2032) to buffer the internal clock in case of power failure.

In normal operation of the nurse call system, i.e. the IP-SystemManager is supplied with power, the battery should be replaced every 2 years. If the IP-SystemManager is not supplied with power, e.g. when the device is stored, an earlier replacement may be necessary due to the battery consumption.

Replacement battery: High-quality Li battery, 3 V, CR2032

**Behaviour of the system when the battery is empty**

An empty battery leads to the following situation: After a failure of the supply voltage of the IP-SystemManager, the internal clock restarts at 00:00. An incorrect time is displayed on the display devices of the nurse call system (ComStation, corridor display, ComTerminals, RoomTerminals). The time is correct again when a connected TimeServer sends the time. If no TimeServer is connected, the time must be reset manually via the SystemOrganizer.

**Battery disposal**

The used battery must be returned to a collection point for used batteries or to Tunstall.

**Technical data**

Dimensions (HxWxD):	132 x 216 x 48 mm
Housing material:	PS UL94-HB
Service cover material:	ABS
Weight:	415 g
IP protection class:	IP 20
Nominal voltage:	24 V DC
Permitted voltage range:	20 – 26 V DC
Ambient temperature:	+5 °C – +40 °C
Relative humidity:	0% – 85% (non condensing)
Standby current consumption:	120 mA

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Câble patch, min. cat. 5e.

## Connections

\*The connections marked with a star are only enabled, if the IP-SystemManager was prepared accordingly at the factory, refer to the project specific documents.

### Connection "Power"

Connection to the 24 V power supply ring of the ward or directly to a 24 V power supply unit. Insert insertion bridges!  
Conductor cross-section: 2.5 mm<sup>2</sup>  
Connector, 4-pole: spare part order no. 00 0211 33.  
Two insertion bridges: spare part order no. 00 0220 52.

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Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>  
Connector, 2-pole: spare part order no. 00 0211 47.

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00 0211 45; Connector, 4-pole, for NF: spare part order no.

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**Green LED "Status"**

LED flashes once a second. = IP-SystemManager is ok.

LED lights up permanently or is off. = Fault! IP-SystemManager is not operational.

**Red LED "Fault" below "Fault Relay"**

LED lights up permanently. = There is at least one fault.

Signalling in parallel to "Fault Relay" output.

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- An OSYnet bus user is defective or disconnected.
- A RAN bus user is defective or disconnected.
- The power supply on a ward is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.
- The power supply at an IP-SystemManager is interrupted This fault is recognized indirectly by a failure in OSYnet bus user polling.
- An OSYnet bus line is interrupted.
- LAN connection to an IP-SystemManager is interrupted.
- An IP-SystemManager is not operational.
- The power supply unit is in buffer mode. Indication in parallel to the red LED "Fault" below "UPS Input".

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© Tunstall GmbH, Orkotten 66, 48291 Teigte, Germany, www.tunstall.de