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**Call/presence switch insert with beeper**

Intended for connection to the room bus of a CONCENTO<sup>CARE</sup> nurse call system.


**NOTE!** The complete installation of the system is described in the Technical Manual.

**A Product description**

The red call button [7] serves for triggering calls or WC calls. An LED in the call button lights up faintly to locate the button in the dark (location light). The LED lights up brightly as soon as a call is triggered (reassurance light).

The green presence button [6] serves for switching the presence of the staff on and off. An LED in the presence button lights up when the presence is switched on (reminder light).

An integrated beeper serves for displaying calls when the presence is switched on. **IMPORTANT!** The beeper is only active if the call/presence switch insert with beeper is connected to a room terminal "control module ZLB" (order no. 19 0701 20). The beeper is inactive with all other room terminals.

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- |  |   |
|--|---|
| [1] *Back box, 1-gang<br>(solid wall: 17 0100 00,<br>partition wall: 17 5100 00) | [5] Terminal block<br>(plug-in terminals) |
| [2] Mounting plate   | [6] Presence button (green)               |
| [3] *Back box screws   | [7] Call button (red)                     |
| [4] *Cover frame, 1-gang   | [8] Switch insert                         |

\* Not included in the scope of delivery, please order separately.

**B Mounting**

Wall mounting on one-gang back box.

1. Bolt the mounting plate [2] to the back box [1]. The four square recesses in the mounting plate are located at the top and bottom.
2. Prepare the wires and feed them through the mounting plate [2] (see section "Connection").
3. Connect the switch insert [8] according to section "Connection".
4. Hold the cover frame [4] in front of the mounting plate [2] on the wall.
5. Plug the switch insert [8] onto the mounting plate so that the latching hook of the switch insert fits into the square recesses of the mounting plate.

**C Dismantling**

1. Pull the cover frame [4] together with the switch insert [8] off the mounting plate [2].
2. For disconnecting the wires, refer to fig. E.

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
**D Connection**

1. Pull the connection cable out through the mounting plate and strip it to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires according to the connection diagram.

**NOTE!** The terminal block [5] can be pulled off temporarily.

3. For the call button [7] to be used for triggering WC calls, you must fit a wire bridge between connection points LT1 and LT2. No wire bridge is to be set for triggering standard calls (= factory setting).

**NOTE!** If the call button is used to trigger WC calls, there must be a facility to cancel the WC call.

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**Technical data**

Nominal voltage	24 V=
Standby current consumption	9 mA
Room bus connection	Cable type: J-Y(St)Y 2x2x0.8
Skinning length	6 mm
Installation depth	24 mm
Degree of protection	IP 22
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 95 %