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Tunstall

Intended use

The connection sockets described in these installation instructions are intended for connection to the RAN room bus with speech in a nurse call system from the Flamenco system family. The connection sockets have a call button as well as a magnetic socket and an 8-pin socket for connecting operating devices.

The connection sockets are designed for wall installation, only in horizontal orientation, or for under-furniture mounting.

Table 1. Products

70042700	Dual connection socket MAG
70042750	Dual connection socket MAG, TVL

Safety notes

NOTICE
The complete installation of the system is described in the technical manual.

CAUTION
The PCB is fitted with electrostatically sensitive components. Therefore, avoid direct contact.

Product description

The **magnetic socket** is used to connect a device with a magnetic plug, e.g.

- ePat® lite MAG (77037010)
- Patient handset MAG (76072000)
- Pear push switch MAG (76071003, 76071103)

The **8-pin socket** is used to connect a device with an 8-pin mini DIN plug, e.g.

- Pear push switch (7007100x, 7007110x, 70071200)


The **red call button** is used to raise calls. The call button lights up dimly to help you find the button in the dark. It lights up brightly as soon as a call is triggered.

The connection supports the following functions of the connected devices: transmission of calls, control of two light sources, TV transmission, Entertainment transmission, blinds control.

NOTICE
When using an ePat® lite, entertainment and blinds control are available. When using a patient handset, either entertainment or blinds control is available.

Figure 1. Control elements

- [1] Call button (red)
- [2] 8-pin socket
- [3] Magnetic socket

Interfaces

Lighting control

WARNING
The use of inadmissible relays can endanger people's lives! Observe the following requirements for the light relays.

The switching outputs LT1 (reading light) and LT2 (room light) provide an output voltage of 24 V DC and a current of max. 40 mA, which is taken from the power supply of the room terminal. The potential refers to the 0 V connection of the room terminal. The switching pulse (24 V DC, 40 mA) supplied by the outputs lasts as long as the button on the control device (ePat® lite, Patient handset, pear push switch) is pressed.

Table 2. Light relay requirements

Application	Switching light on/off (without dimming)	Dimming the light
Relay type	• Pulse relay (electronic) • Pulse relay (mechanical)	• Switching relay (electronic) • Switching relay (mechanical)
Nominal control voltage	24 V DC	
Control voltage range	18 – 26 V DC	
Maximum current consumption	40 mA	
Freewheeling diode	 When using mechanical relays, a freewheeling diode (e.g. 1N4002) directly on the relay is essential.	
Potential separation	Safe isolation with 2 x MOPP according to DIN EN 60601-1 must be provided.	
Ballast	-	For the dimming function, a corresponding dimmable ballast is required (e.g. OS-RAM DALI).

TV control

The television set is not connected to the Dual connection socket MAG, but via the RAN room bus, e.g. via a RAN interface (article no. 77084000) or via a IR TV control module universal (article no. 77036011).


Blinds control

The blinds control is not connected to the Dual connection socket MAG, but to a RAN interface Universal (article no. 70084800), which is connected to the RAN room bus.

The relays connected to the RAN interface Universal must be self-powered, controllable by a normally open contact and have a safe separation with 2 x MOPP in accordance with DIN EN 60601-1. Maximum contact load: 60 mA / 24 V.

Setting DIP switches

Figure 2. DIP switch on the PCB

NOTE! DIP switch **PG** (4) must always be set to OFF.

Setting the bed number

A bed number must be set on each connection socket:

Bed number	DIP switch		
	P3 (1)	P2 (2)	P1 (3)
Bed 1	OFF	OFF	ON
Bed 2	OFF	ON	OFF
Bed 3	OFF	ON	ON
Bed 4	ON	OFF	OFF
Bed 5	ON	OFF	ON
Bed 6	ON	ON	OFF


Setting blinds option

Blinds option	TV (5) DIP switch
No blinds control connected.	OFF
Blinds control	ON

Mounting

The connection sockets are designed for wall installation, only in horizontal orientation, or for under-furniture mounting.

Figure 3. Mounting

- [1] Socket element
- [2] Mounting frame
- [3] Two plugs of flat cables
- [4] Cover plate with PCB
- [5] Socket for plug A
- [6] Socket for plug B
- [7] DIP switch (5-pole)


1. Insert the socket element into the back box.

2. Screw the mounting frame together with the socket element onto the back box.
3. Insert plug **A** and plug **B** (see Figure 6) into the corresponding sockets.
4. Press the cover plate onto the mounting frame.

Dismantling

1. Lever the cover plate with built-in PCB off the wall with both hands until it comes loose.
2. Disconnect both flat cable plugs from the PCB so that you can remove the cover plate with the built-in PCB completely from the wall.

Figure 4. Dismantling

Technical data


Dimensions (HxWxD)	90 x 180 x 45 mm
Mounting depth	30 mm
Housing material	Polycarbonate
Housing colour	Studio white, similar to RAL 9016
Power supply	24 V DC
Quiescent current consumption	39 mA
Max. current consumption (without connected devices)	51 mA
TV sound input from 70042700	25 Veff
TV sound input from 70042750	5 Veff
Ambient temperature	+5 °C to +40 °C
Relative humidity	0 % to 85 % (non-condensing)

Tunstall GmbH, Orkotten 66, 48291 Telgte, Germany
 Tel: +49 (0) 2504/701-0, E-Mail: DE.Info@tunstall.com
 www.tunstall.de/en-de/

Tunstall


Connections

Figure 5. Connections on the connection block

Notes on the connections

Figure 6. Connection block

Handling the push wire connectors

Usable conductor cross sections: 0.55 - 1 mm Ø.

Connect only one wire per opening.


1. Strip the wire (8 mm).
2. Insert the stripped wire into the terminal until it stops.
3. To disconnect the conductor, insert a screwdriver (blade: 2.5 x 0.4 mm or 3.5 x 0.5 mm) into the opening and pull out the conductor.

Connect shield wire

To connect a shield wire, insert a screwdriver (blade: 2.5 x 0.4 mm or 3.5 x 0.5 mm) into the opening and hold the spring contact open with the screwdriver while inserting the shield wire.

Then check that the contact is correctly connected by pulling on the conductor.

Figure 7. Connect shield wire

RAN room bus

- The cable for speech transmission, i.e. the cable containing the wires **MIC**, **LS**, and **SCH-B**, may be stripped for a max. of 30 mm only, to ensure the protective screen function. Otherwise, a short circuit situation or crosstalk may occur.
- The cable shield and the shield wire of the cable containing the wires **RAN**, **+24V-Si** and **0V** are not connected and must be removed in the back box up to the cable sheath.
- The unused wires must not be interconnected.

Tunstall GmbH, Orkotten 66, 48291 Telgte, Germany
Tel: +49 (0) 2504/701-0, E-Mail: DE.Info@tunstall.com
www.tunstall.de/en-de/

Tunstall